

Beta-defensin 119

Catalog No: #AP79719



Package Size: #AP79719-1 50ug #AP79719-2 100ug #AP79719-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-defensin 119
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Beta-defensin 120,Beta-defensin 19,DEFB-19,Beta-defensin 20,DEFB-20,Defensin, beta 119,Defensin, beta 120,ESC42-RELA,DEFB120,DEFB19,DEFB20,UNQ2449/PRO5729
Accession No.	Uniprot ID:Q8N690
Uniprot	Q8N690
GeneID	245932;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DEFB119

Note: This product is for in vitro research use only