

GST Monoclonal Antibody(Q16)

Catalog No: #44023

Package Size: #44023-1 50ul #44023-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GST Monoclonal Antibody(Q16)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Q16
Isotype	IgG1
Applications	WB
Specificity	The antibody detects C-terminal, internal, and N-terminal GST fusion proteins.
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	GST
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:2000

Images



Western blot analysis of 0.2ug GST protein, diluted at 1:5000

Note: This product is for in vitro research use only