

Matrix metalloproteinase-26

Catalog No: #AP79722



Package Size: #AP79722-1 50ug #AP79722-2 100ug #AP79722-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Matrix metalloproteinase-26
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MMP-26,Endometase,Matrilysin-2
Accession No.	Uniprot ID:Q9NRE1
Uniprot	Q9NRE1
GeneID	56547;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MMP26

Note: This product is for in vitro research use only