

Cocaine esterase

Catalog No: #AP79723



Package Size: #AP79723-1 50ug #AP79723-2 100ug #AP79723-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cocaine esterase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Carboxylesterase 2,CE-2,Methylumbelliferyl-acetate deacetylase 2,ICE
Accession No.	Uniprot ID:O00748
Uniprot	O00748
GeneID	8824;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CES2

Note: This product is for in vitro research use only