

# Leucine-rich repeat-containing G-protein coupled receptor 4

Catalog No: #AP79726

Package Size: #AP79726-1 50ug #AP79726-2 100ug #AP79726-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Leucine-rich repeat-containing G-protein coupled receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	G-protein coupled receptor 48,GPR48
Accession No.	Uniprot ID:Q9BXB1
Uniprot	Q9BXB1
GeneID	55366;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:LGR4

Note: This product is for in vitro research use only