

Baculoviral IAP repeat-containing protein 5

Catalog No: #AP79727



Package Size: #AP79727-1 50ug #AP79727-2 100ug #AP79727-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Baculoviral IAP repeat-containing protein 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Apoptosis inhibitor 4, Apoptosis inhibitor survivin, API4, IAP4
Accession No.	Uniprot ID:O15392
Uniprot	O15392
GeneID	332;
SDS-PAGE MW	15.62kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:1-142aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BIRC5

Note: This product is for in vitro research use only