

Bifunctional heparan sulfate N-deacetylase/N-sulfotransferase 2

Catalog No: #AP79729

Package Size: #AP79729-1 50ug #AP79729-2 100ug #AP79729-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bifunctional heparan sulfate N-deacetylase/N-sulfotransferase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Glucosaminyl N-deacetylase/N-sulfotransferase 2,NDST-2,N-heparan sulfate sulfotransferase 2,N-HSST 2,HSST2
Accession No.	Uniprot ID:P52849
Uniprot	P52849
GeneID	8509;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NDST2

Note: This product is for in vitro research use only