

26S protease regulatory subunit 6A

Catalog No: #AP79730



Package Size: #AP79730-1 50ug #AP79730-2 100ug #AP79730-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	26S protease regulatory subunit 6A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	26S proteasome AAA-ATPase subunit RPT5, Proteasome 26S subunit ATPase 3, Proteasome subunit P50, Tat-binding protein 1, TBP-1, TBP1
Accession No.	Uniprot ID:P17980
Uniprot	P17980
GeneID	5702;
SDS-PAGE MW	16.72kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:1-152aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PSMC3

Note: This product is for in vitro research use only