

# Inactive tyrosine-protein kinase transmembrane receptor ROR1

Catalog No: #AP79732

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP79732-1 50ug #AP79732-2 100ug #AP79732-3 1mg

## Description

Product Name	Inactive tyrosine-protein kinase transmembrane receptor ROR1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Neurotrophic tyrosine kinase, receptor-related 1,NTRKR1
Accession No.	Uniprot ID:Q01973
Uniprot	Q01973
GeneID	4919;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ROR1

Note: This product is for in vitro research use only