

# Paired immunoglobulin-like type 2 receptor beta

Catalog No: #AP79733



Package Size: #AP79733-1 50ug #AP79733-2 100ug #AP79733-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Paired immunoglobulin-like type 2 receptor beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Activating receptor PILR-beta, Cell surface receptor FDFACT, FDFACT, PP1551
Accession No.	Uniprot ID:Q9UKJ0
Uniprot	Q9UKJ0
GeneID	29990;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: PILRB

Note: This product is for in vitro research use only