

TNF receptor-associated factor 4

Catalog No: #AP79736



Package Size: #AP79736-1 50ug #AP79736-2 100ug #AP79736-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TNF receptor-associated factor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cysteine-rich domain associated with RING and Traf domains protein 1, Metastatic lymph node gene 62 protein, MLN 62, RING finger protein 83, CART1, MLN62, RNF83
Accession No.	Uniprot ID: Q9BUZ4
Uniprot	Q9BUZ4
GeneID	9618;
SDS-PAGE MW	21.34kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region: 1-194aa Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TRAF4

Note: This product is for in vitro research use only