

Myosin-binding protein C, cardiac-type

Catalog No: #AP79737



Package Size: #AP79737-1 50ug #AP79737-2 100ug #AP79737-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myosin-binding protein C, cardiac-type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Cardiac MyBP-C,C-protein, cardiac muscle isoform
Accession No.	Uniprot ID:Q90688
Uniprot	Q90688
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MYBPC3

Note: This product is for in vitro research use only