

Heat shock protein beta-1

Catalog No: #AP79738



Package Size: #AP79738-1 50ug #AP79738-2 100ug #AP79738-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heat shock protein beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HspB1,28 kDa heat shock protein,Estrogen-regulated 24 kDa protein,Heat shock 27 kDa protein,HSP 27,Stress-responsive protein 27,SRP27,HSP27,HSP28
Accession No.	Uniprot ID:P04792
Uniprot	P04792
GeneID	3315;
SDS-PAGE MW	22.55kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:1-205aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HSPB1

Note: This product is for in vitro research use only