

Cytochrome P450 2D18

Catalog No: #AP79739



Package Size: #AP79739-1 50ug #AP79739-2 100ug #AP79739-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome P450 2D18
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CYP11D18,Cytochrome P450 2D-29,Cytochrome P450 2D-35,Cyp2d-18
Accession No.	Uniprot ID:Q64680
Uniprot	Q64680
GeneID	171522;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cyp2d18

Note: This product is for in vitro research use only