

B-lymphocyte antigen CD20

Catalog No: #AP79741



Package Size: #AP79741-1 50ug #AP79741-2 100ug #AP79741-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	B-lymphocyte antigen CD20
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	B-lymphocyte surface antigen B1,Bp35,Leukocyte surface antigen Leu-16,Membrane-spanning 4-domains subfamily A member 1,CD20,CD20
Accession No.	Uniprot ID:P11836
Uniprot	P11836
GeneID	931;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MS4A1

Note: This product is for in vitro research use only