

# Proliferating cell nuclear antigen

Catalog No: #AP79742



Package Size: #AP79742-1 50ug #AP79742-2 100ug #AP79742-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Proliferating cell nuclear antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PCNA,Cyclin
Accession No.	Uniprot ID:P12004
Uniprot	P12004
GeneID	5111;
SDS-PAGE MW	28.71kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Expression Region:1-261aaStorage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PCNA

Note: This product is for in vitro research use only