

Gamma-interferon-inducible lysosomal thiol reductase

Catalog No: #AP79744



Package Size: #AP79744-1 50ug #AP79744-2 100ug #AP79744-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gamma-interferon-inducible lysosomal thiol reductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Gamma-interferon-inducible protein IP-30,Legumaturain,GILT,IP30
Accession No.	Uniprot ID:P13284
Uniprot	P13284
GeneID	10437;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IFI30

Note: This product is for in vitro research use only