

Phospholysine phosphohistidine inorganic pyrophosphate phosphatase

Catalog No: #AP79745

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP79745-1 50ug #AP79745-2 100ug #AP79745-3 1mg

Description

Product Name	Phospholysine phosphohistidine inorganic pyrophosphate phosphatase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	hLHPP
Accession No.	Uniprot ID:Q9H008
Uniprot	Q9H008
GeneID	64077;
SDS-PAGE MW	29.7kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:1-270aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LHPP

Note: This product is for in vitro research use only