

Delta-like protein 1

Catalog No: #AP79746



Package Size: #AP79746-1 50ug #AP79746-2 100ug #AP79746-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Delta-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Drosophila Delta homolog 1,Delta1,UNQ146/PRO172
Accession No.	Uniprot ID:O00548
Uniprot	O00548
GeneID	28514;
SDS-PAGE MW	23.65kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:131-345aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DLL1

Note: This product is for in vitro research use only