

D-dopachrome decarboxylase

Catalog No: #AP79747



Package Size: #AP79747-1 50ug #AP79747-2 100ug #AP79747-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	D-dopachrome decarboxylase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	D-dopachrome tautomerase
Accession No.	Uniprot ID:A5PK65
Uniprot	A5PK65
GeneID	615999;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: DDT

Note: This product is for in vitro research use only