

## Phosphoglycerate kinase 1

Catalog No: #AP79748



Package Size: #AP79748-1 50ug #AP79748-2 100ug #AP79748-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phosphoglycerate kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cell migration-inducing gene 10 protein,Primer recognition protein 2,PRP 2,PGKA,MIG10,OK/SW-cl.110
Accession No.	Uniprot ID:P00558
Uniprot	P00558
GeneID	5230;
SDS-PAGE MW	45.87kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Expression Region:1-417aaStorage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PGK1

Note: This product is for in vitro research use only