

S-formylglutathione hydrolase

Catalog No: #AP79750



Package Size: #AP79750-1 50ug #AP79750-2 100ug #AP79750-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	S-formylglutathione hydrolase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FGH,Esterase 10,Esterase D,Sid 478,Es10,Sid478
Accession No.	Uniprot ID:Q9R0P3
Uniprot	Q9R0P3
GeneID	13885;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Esd

Note: This product is for in vitro research use only