

Vitronectin

Catalog No: #AP79752



Package Size: #AP79752-1 50ug #AP79752-2 100ug #AP79752-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vitronectin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	VN,S-protein,Serum-spreading factor
Accession No.	Uniprot ID:P29788
Uniprot	P29788
GeneID	22370;
SDS-PAGE MW	36.08kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:151-478aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Vtn

Note: This product is for in vitro research use only