

Methylcytosine dioxygenase TET3

Catalog No: #AP79754



Package Size: #AP79754-1 50ug #AP79754-2 100ug #AP79754-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Methylcytosine dioxygenase TET3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:Q8BG87
Uniprot	Q8BG87
GeneID	194388;
SDS-PAGE MW	12.54kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:1555-1668aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tet3

Note: This product is for in vitro research use only