

Pulmonary surfactant-associated protein B

Catalog No: #AP79756



Package Size: #AP79756-1 50ug #AP79756-2 100ug #AP79756-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pulmonary surfactant-associated protein B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	SP-B,6 kDa protein,Pulmonary surfactant-associated proteolipid SPL(Phe),SFTP3
Accession No.	Uniprot ID:P15781
Uniprot	P15781
GeneID	507398;
SDS-PAGE MW	17.6kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:177-336aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SFTPB

Note: This product is for in vitro research use only