

Prominin-1

Catalog No: #AP79757



Package Size: #AP79757-1 50ug #AP79757-2 100ug #AP79757-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prominin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Antigen AC133,Prominin-like protein 1,CD133,PROML1,MSTP061
Accession No.	Uniprot ID:O43490
Uniprot	O43490
GeneID	8842;
SDS-PAGE MW	28.05kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:179-433aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PROM1

Note: This product is for in vitro research use only