

Islet amyloid polypeptide

Catalog No: #AP79760



Package Size: #AP79760-1 50ug #AP79760-2 100ug #AP79760-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Islet amyloid polypeptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Amylin,Diabetes-associated peptide,DAP,Insulinoma amyloid peptide
Accession No.	Uniprot ID:P10997
Uniprot	P10997
GeneID	3375;
SDS-PAGE MW	9.79kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:1-89aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IAPP

Note: This product is for in vitro research use only