

LIM domain and actin-binding protein 1

Catalog No: #AP79761



Package Size: #AP79761-1 50ug #AP79761-2 100ug #AP79761-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LIM domain and actin-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Epithelial protein lost in neoplasm,mEPLIN,D15Erd366e,Eplin
Accession No.	Uniprot ID:Q9ERG0
Uniprot	Q9ERG0
GeneID	65970;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lima1

Note: This product is for in vitro research use only