

Disintegrin and metalloproteinase domain-containing protein 12

Catalog No: #AP79762

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP79762-1 50ug #AP79762-2 100ug #AP79762-3 1mg

Description

| | |
|--------------------|--|
| Product Name | Disintegrin and metalloproteinase domain-containing protein 12 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | ADAM 12,Meltrin-alpha,Mltna |
| Accession No. | Uniprot ID:Q61824 |
| Uniprot | Q61824 |
| GeneID | 11489; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Adam12

Note: This product is for in vitro research use only