

# Programmed cell death 1 ligand 1

Catalog No: #AP79764

Package Size: #AP79764-1 50ug #AP79764-2 100ug #AP79764-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Programmed cell death 1 ligand 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PD-L1,B7 homolog 1,B7-H1,CD274,B7H1,PDCD1L1,PDCD1LG1,PDL1
Accession No.	Uniprot ID:Q9NZQ7
Uniprot	Q9NZQ7
GeneID	29126;
SDS-PAGE MW	29.92kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Expression Region:19-290aaStorage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CD274

Note: This product is for in vitro research use only