

RuvB-like 2

Catalog No: #AP79770



Package Size: #AP79770-1 50ug #AP79770-2 100ug #AP79770-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RuvB-like 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	48 kDa TATA box-binding protein-interacting protein, 48 kDa TBP-interacting protein, 51 kDa erythrocyte cytosolic protein, ECP-51, INO80 complex subunit J, Repressing pontin 52, Reptin 52, TIP49b, TIP60-associated protein 54-beta, TAP54-beta, INO80J, TIP48, TIP49B, CGI-46
Accession No.	Uniprot ID: Q9Y230
Uniprot	Q9Y230
GeneID	10856;
SDS-PAGE MW	19.25kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region: 210-384aa Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RUVBL2

Note: This product is for in vitro research use only