

Interleukin-2

Catalog No: #AP79775



Package Size: #AP79775-1 50ug #AP79775-2 100ug #AP79775-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	IL-2,T-cell growth factor,TCGF,IL-2
Accession No.	Uniprot ID:P05016
Uniprot	P05016
GeneID	280822;
SDS-PAGE MW	14.85kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:21-155aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL2

Note: This product is for in vitro research use only