

Non-receptor tyrosine-protein kinase TNK1

Catalog No: #AP79777



Package Size: #AP79777-1 50ug #AP79777-2 100ug #AP79777-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Non-receptor tyrosine-protein kinase TNK1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CD38 negative kinase 1
Accession No.	Uniprot ID:Q13470
Uniprot	Q13470
GeneID	8711;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TNK1

Note: This product is for in vitro research use only