

Lysophospholipid acyltransferase 5

Catalog No: #AP79778



Package Size: #AP79778-1 50ug #AP79778-2 100ug #AP79778-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysophospholipid acyltransferase 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LPLAT 5,1-acylglycerophosphocholine O-acyltransferase,1-acylglycerophosphoethanolamine O-acyltransferase,1-acylglycerophosphoserine O-acyltransferase,Lysophosphatidylcholine acyltransferase,LPCAT,Lysophosphatidylcholine acyltransferase 3,Lyso-PC acyltransferase 3,Lysophosphatidylethanolamine acyltransferase,LPEAT,Lysophosphatidylserine acyltransferase,LPSAT,Membrane-bound O-acyltransferase domain-containing protein 5,O-acyltransferase domain-containing protein 5,Grc3f,Mboat5,Oact5
Accession No.	Uniprot ID:Q91V01
Uniprot	Q91V01
GeneID	14792;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lpcat3

Note: This product is for in vitro research use only