

## Angiotensin-converting enzyme 2

Catalog No: #AP79779

Package Size: #AP79779-1 50ug #AP79779-2 100ug #AP79779-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Angiotensin-converting enzyme 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Angiotensin-converting enzyme 2,Angiotensin-converting enzyme 2,Angiotensin-converting enzyme 2
Accession No.	Uniprot ID:Q8NI99
Uniprot	Q8NI99
GeneID	83854;
SDS-PAGE MW	24.31kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Expression Region:221-441aaStorage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ANGPTL6

Note: This product is for in vitro research use only