

Aquaporin-9

Catalog No: #AP79780



Package Size: #AP79780-1 50ug #AP79780-2 100ug #AP79780-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Aquaporin-9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AQP-9,Aquaglyceroporin-9,Small solute channel 1,SSC1
Accession No.	Uniprot ID:O43315
Uniprot	O43315
GeneID	366;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AQP9

Note: This product is for in vitro research use only