

Chemokine-like receptor 1

Catalog No: #AP79782



Package Size: #AP79782-1 50ug #AP79782-2 100ug #AP79782-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chemokine-like receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	G-protein coupled receptor ChemR23,G-protein coupled receptor DEZ,CHEMR23,DEZ
Accession No.	Uniprot ID:Q99788
Uniprot	Q99788
GeneID	1240;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CMKLR1

Note: This product is for in vitro research use only