

Gremlin-1

Catalog No: #AP79786



Package Size: #AP79786-1 50ug #AP79786-2 100ug #AP79786-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Gremlin-1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Cysteine knot superfamily 1, BMP antagonist 1,Down-regulated in Mos-transformed cells protein,Cktsf1b1,Drm |
| Accession No. | Uniprot ID:O70326 |
| Uniprot | O70326 |
| GeneID | 23892; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Grem1

Note: This product is for in vitro research use only