

Semaphorin-3E

Catalog No: #AP79788



Package Size: #AP79788-1 50ug #AP79788-2 100ug #AP79788-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Semaphorin-3E
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Semaphorin-H,Sema H,Semah,Semh
Accession No.	Uniprot ID:P70275
Uniprot	P70275
GeneID	20349;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sema3e

Note: This product is for in vitro research use only