

Filamin-B

Catalog No: #AP79790



Package Size: #AP79790-1 50ug #AP79790-2 100ug #AP79790-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Filamin-B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FLN-B,ABP-280-like protein,Actin-binding-like protein,Beta-filamin
Accession No.	Uniprot ID:Q80X90
Uniprot	Q80X90
GeneID	286940;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Flnb

Note: This product is for in vitro research use only