

Fibromodulin

Catalog No: #AP79791



Package Size: #AP79791-1 50ug #AP79791-2 100ug #AP79791-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Fibromodulin |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | FM, Collagen-binding 59 kDa protein, Keratan sulfate proteoglycan fibromodulin, KSPG fibromodulin, FM, SLRR2E |
| Accession No. | Uniprot ID:Q06828 |
| Uniprot | Q06828 |
| GeneID | 2331; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FMOD

Note: This product is for in vitro research use only