

# 6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 4

Catalog No: #AP79792

Package Size: #AP79792-1 50ug #AP79792-2 100ug #AP79792-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	6PF-2-K/Fru-2,6-P2ase 4,6PF-2-K/Fru-2,6-P2ase testis-type isozyme
Accession No.	Uniprot ID:P25114
Uniprot	P25114
GeneID	100911725;54283;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pfkfb4

Note: This product is for in vitro research use only