

Collagen alpha-1(II) chain

Catalog No: #AP79794



Package Size: #AP79794-1 50ug #AP79794-2 100ug #AP79794-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Collagen alpha-1(II) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Alpha-1 type II collagen
Accession No.	Uniprot ID:AAH07252.1/ncbi
Uniprot	AAH07252.1
SDS-PAGE MW	27.06kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:23-268aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:COL2A1

Note: This product is for in vitro research use only