

Fibrinogen-like protein 1

Catalog No: #AP79795



Package Size: #AP79795-1 50ug #AP79795-2 100ug #AP79795-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibrinogen-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Mouse fibrinogen-related protein 1,MFIRE-1,Mfire1
Accession No.	Uniprot ID:Q71KU9
Uniprot	Q71KU9
GeneID	234199;
SDS-PAGE MW	32.12kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:23-314aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fgl1

Note: This product is for in vitro research use only