

## Enhancer of filamentation 1

Catalog No: #AP79796



Package Size: #AP79796-1 50ug #AP79796-2 100ug #AP79796-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Enhancer of filamentation 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	hEF1,CRK-associated substrate-related protein,CAS-L,Cas scaffolding protein family member 2,Neural precursor cell expressed developmentally down-regulated protein 9,NEDD-9,Renal carcinoma antigen NY-REN-12,p105,CASL,CASS2
Accession No.	Uniprot ID:Q14511
Uniprot	Q14511
GeneID	4739;
SDS-PAGE MW	20.02kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Expression Region:234-415aaStorage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:NEDD9

Note: This product is for in vitro research use only