

## Anionic trypsin-2

Catalog No: #AP79802

Package Size: #AP79802-1 50ug #AP79802-2 100ug #AP79802-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Anionic trypsin-2   |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE                                |
| Species Reactivity | Rt  |
| Other Names        | Anionic trypsin II,Pretrypsinogen II,Serine protease 2,Try2 |
| Accession No.      | Uniprot ID:P00763   |
| Uniprot            | P00763  |
| GeneID             | 25052;  |
| SDS-PAGE MW        | 24.53kDa  |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)      |

## Application Details

Expression Region:24-246aaStorage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Prss2

Note: This product is for in vitro research use only