

Non-receptor tyrosine-protein kinase TYK2

Catalog No: #AP79803



Package Size: #AP79803-1 50ug #AP79803-2 100ug #AP79803-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Non-receptor tyrosine-protein kinase TYK2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:P29597
Uniprot	P29597
GeneID	7297;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TYK2

Note: This product is for in vitro research use only