

## Sugar transporter SWEET1

Catalog No: #AP79804



Package Size: #AP79804-1 50ug #AP79804-2 100ug #AP79804-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Sugar transporter SWEET1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HsSWEET1,RAG1-activating protein 1,Solute carrier family 50 member 1,Stromal cell protein,RAG1AP1,SCP
Accession No.	Uniprot ID:Q9BRV3
Uniprot	Q9BRV3
GeneID	55974;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SLC50A1

Note: This product is for in vitro research use only