

## CAAX prenyl protease 1 homolog

Catalog No: #AP79806



Package Size: #AP79806-1 50ug #AP79806-2 100ug #AP79806-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CAAX prenyl protease 1 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Farnesylated proteins-converting enzyme 1,FACE-1,Prenyl protein-specific endoprotease 1,Zinc metalloproteinase Ste24 homolog,FACE1,STE24
Accession No.	Uniprot ID:O75844
Uniprot	O75844
GeneID	10269;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ZMPSTE24

Note: This product is for in vitro research use only