

# Kynurenine--oxoglutarate transaminase 1, mitochondrial

Catalog No: #AP79808

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP79808-1 50ug #AP79808-2 100ug #AP79808-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Kynurenine--oxoglutarate transaminase 1, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Kynurenine--oxoglutarate transaminase I, Cysteine-S-conjugate beta-lyase, Glutamine transaminase K, GTK, Glutamine--phenylpyruvate transaminase, Kynurenine aminotransferase 1, Kynurenine aminotransferase I, KATI, Ccb11, Kat
Accession No.	Uniprot ID:Q08415
Uniprot	Q08415
GeneID	311844;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Kyat1

Note: This product is for in vitro research use only